



Developed By P-I Brånemark

Instrument Guidance



A COMPANY OF
SIC invent Group

P-I Instrument Kit System Reference Guide

General Cleaning Instructions

GENERAL INSTRUCTIONS

- Instruments must be cleaned and sterilized individually, without coming into contact with any other instrument
- Always inspect Instruments before use, assure there are no mechanical damages and tools are properly assembled
- Surgical Instruments supplied as non sterile are to be sterilized before use
- Clean and rinse your Instrument before sterilizing
- Do not try to regrind or refurbish
- Dispose the device after its lifetime in accordance to local sanitary legislation
- Instruments are to be used by qualified and trained professionals only

CLEANING INSTRUCTIONS

Cleaning of Instruments

- Clean manually using mild soap solution with a soft bristle toothbrush to remove any stains. Ensure no traces of detergents adhere to the Instruments before dipping into the disinfectant solution. Use disinfectant solutions according to the manufacturer's instructions
- Do not allow saline solutions or any protein residue to dry on Instruments
- Alternatively, ultrasonic cleaning method can be used. If used, strictly adhere to the ultrasonic manufacturer's instructions while using this method. Use specified cleaning solutions. Change the solutions frequently or as instructed by the manufacturer

Cleaning of Trays | Cases

NOTE: the information below is applicable to polymer (RADEL R polyphenylsulfone resin) and stainless steel surgical trays.

- Remove all parts and inserts from the surgical tray
- Rinse the tray and tray insert with cool to lukewarm tap water
- Use a damp cloth to wipe and remove any excess soil from each part
- Wipe each part with a cloth dipped in an enzymatic detergent diluted with tap water per the manufacturer's instructions
- Rinse thoroughly with tap water for 3 min

STERILIZATION

NOTE: Chemclave sterilization is NOT recommended.

- Dry parts thoroughly before sterilization
- Double-wrap kit with polypropylene wrap and secure
- wrap with autoclave tape
- Recommended Sterilization Parameters:

CYCLE TYPE (STEAM)	TEMPERATURE	EXPOSURE TIME	DRY TIME (MIN.)
PRE-VACUUM ^{1, 2}	132°C/270°F	4 minutes	20 minutes
GRAVITY ¹	121°C	55 minutes	20 minutes
GRAVITY	132°C	55 minutes	20 minutes

¹ Minimum validated sterilization time and temperature required to achieve a 10-6 sterility assurance level (SAL).
² Local or national specifications should be followed where steam sterilization requirements are stricter or more conservative than those listed in this table.

Instrument Kit System Reference Guide



POLYMER KIT

Advanced

Application:	Surgical
	Prosthetic
	MT (ALL)
	AMP-S and AMP (ALL)
	HEX-S and HEX (ALL)

Content:

Tray	Code
Polymer	131134

Implant Drivers

Implant Insertion Driver - Medium	131139
Implant Insertion Driver - Long	131140
Implant Insertion Driver - HEX 3.5 Medium	131141
Implant Insertion Driver - HEX 3.5 - Long	131142

Drills

First Drill + 2.2	141138
Conical Drill - 2.8	141146
Conical Drill - 3.4	141148
Conical Drill - 3.8	141314
Conical Drill - 4.6	141152
Conical Drill - 4.8	141315
Dense Drill - 3.3	141213
Dense Drill - 3.75	141316
Dense Drill - 4.0	141215
Dense Drill - 4.8 5.0	141317

Instruments

Torque Wrench - Surgical and Prosthetic	131128
Driver Adapter - Prosthetic Squared	131129
Wrench Adapter - Internal Grip	131130
Hexagonal Driver - Medium - 1.2	131011
Conical Abutment Driver - Medium	131017
Retriever - MT	131131

Accessories

Drill Extension	131028
Guide Pin 2.2 2.8 - Conical (2)	131114
Guide Pin 2.8 3.8 - Conical (2)	131115

Start-up

Application:	Surgical
	Prosthetic
	MT Ø3.75
	AMP-S and AMP Ø3.75 and 4.0
	HEX-S and HEX Ø3.75 and 4.0

Content:

Tray	Code
Polymer	131134

Implant Drivers

Implant Insertion Driver - Medium	131139
Implant Insertion Driver - Long	131140

Drills

First Drill + 2.2	141138
Conical Drill - 2.8	141146
Conical Drill - 3.4	141148
Conical Drill - 3.8	141314
Dense Drill - 3.75	141316
Dense Drill - 4.0	141215

Instruments

Torque Wrench - Surgical and Prosthetic	131128
Driver Adapter - Prosthetic Squared	131129
Wrench Adapter - Internal Grip	131130
Hexagonal Driver - Medium - 1.2	131011
Conical Abutment Driver - Medium	131017
Retriever - MT	131131

Accessories

Drill Extension	131028
Guide Pin 2.2 2.8 - Conical (2)	131114
Guide Pin 2.8 3.8 - Conical (2)	131115

POLYMER KIT

Code 131134





STAINLESS STEEL

Advanced

Application:	Surgical
	Prosthetic
	MT (ALL)
	AMP-S and AMP (ALL)
	HEX-S and HEX (ALL)

Content:

Tray	Code
Stainless Steel	131117

Implant Drivers

Implant Insertion Driver - Medium	131139
Implant Insertion Driver - Long	131140
Implant Insertion Driver - HEX 3.5 Medium	131141
Implant Insertion Driver - HEX 3.5 - Long	131142

Drills

First Drill + 2.2	141138
Conical Drill - 2.8	141146
Conical Drill - 3.4	141148
Conical Drill - 3.8	141314
Conical Drill - 4.6	141152
Conical Drill - 4.8	141315
Dense Drill - 3.3	141213
Dense Drill - 3.75	141316
Dense Drill - 4.0	141215
Dense Drill - 4.8 5.0	141317

Instruments

Torque Wrench - Surgical and Prosthetic	131128
Driver Adapter - Prosthetic Squared	131129
Wrench Adapter - Internal Grip	131130
Hexagonal Driver - Medium - 1.2	131011
Conical Abutment Driver - Medium	131017
Retriever - MT	131131

Accessories

Drill Extension	131028
Guide Pin 2.2 2.8 - Conical (2)	131114
Guide Pin 2.8 3.8 - Conical (2)	131115

Start-up

Application:	Surgical
	Prosthetic
	MT Ø3.75
	AMP-S and AMP Ø3.75 and 4.0
	HEX-S and HEX Ø3.75 and 4.0

Content:

Tray	Code
Stainless Steel	131117

Implant Drivers

Implant Insertion Driver - Medium	131139
Implant Insertion Driver - Long	131140

Drills

First Drill + 2.2	141138
Conical Drill - 2.8	141146
Conical Drill - 3.4	141148
Conical Drill - 3.8	141314
Dense Drill - 3.75	141316
Dense Drill - 4.0	141215

Instruments


Torque Wrench - Surgical and Prosthetic	131128
Driver Adapter - Prosthetic Squared	131129
Wrench Adapter - Internal Grip	131130
Hexagonal Driver - Medium - 1.2	131011
Conical Abutment Driver - Medium	131017
Retriever - MT	131131

Accessories

Drill Extension	131028
Guide Pin 2.2 2.8 - Conical (2)	131114
Guide Pin 2.8 3.8 - Conical (2)	131115

STAINLESS STEEL

Code 131117





COMPACT

Advanced

Application:	Surgical
	Prosthetic
	MT (ALL)
	AMP-S and AMP (ALL)
	HEX-S and HEX (ALL)

Content:

Tray	Code
Compact	131138

Implant Drivers

Implant Insertion Driver - Medium	131139
Implant Insertion Driver - Long	131140
Implant Insertion Driver - HEX 3.5 Medium	131141
Implant Insertion Driver - HEX 3.5 - Long	131142

Drills

First Drill + 2.2	141138
Conical Drill - 2.8	141146
Conical Drill - 3.4	141148
Conical Drill - 3.8	141314
Conical Drill - 4.6	141152
Conical Drill - 4.8	141315
Dense Drill - 3.3	141213
Dense Drill - 3.75	141316
Dense Drill - 4.0	141215
Dense Drill - 4.8 5.0	141317

Instruments

Torque Wrench - Surgical and Prosthetic	131128
Driver Adapter - Prosthetic Squared	131129
Wrench Adapter - Internal Grip	131130
Hexagonal Driver - Medium - 1.2	131011
Conical Abutment Driver - Medium	131017
Retriever - MT	131131

Accessories

Drill Extension	131028
Guide Pin 2.2 2.8 - Conical (2)	131114
Guide Pin 2.8 3.8 - Conical (2)	131115

Start-up

Application:	Surgical
	Prosthetic
	MT Ø3.75
	AMP-S and AMP Ø3.75 and 4.0
	HEX-S and HEX Ø3.75 and 4.0

Content:

Tray	Code
Compact	131138

Implant Drivers

Implant Insertion Driver - Medium	131139
Implant Insertion Driver - Long	131140

Drills

First Drill + 2.2	141138
Conical Drill - 2.8	141146
Conical Drill - 3.4	141148
Conical Drill - 3.8	141314
Dense Drill - 3.75	141316
Dense Drill - 4.0	141215

Instruments

Torque Wrench - Surgical and Prosthetic	131128
Driver Adapter - Prosthetic Squared	131129
Wrench Adapter - Internal Grip	131130
Hexagonal Driver - Medium - 1.2	131011
Conical Abutment Driver - Medium	131017
Retriever - MT	131131

Accessories

Drill Extension	131028
Guide Pin 2.2 2.8 - Conical (2)	131114
Guide Pin 2.8 3.8 - Conical (2)	131115



SPECIALIST KITS

HEX | HEX-S

Application: Surgical
HEX-S and HEX
Ø3.75 or 4.0 (a)

Content:	
Tray	Code
Compact	131138

Implant Drivers	
Implant Insertion Driver - Medium	131139
Implant Insertion Driver - Medium-Manual	131146

Drills	
First Drill + 2.2	141138
Conical Drill - 2.8	141146
Conical Drill - 3.4 (a)	141148
Conical Drill - 3.8 (a)	141314
Dense Drill - 3.75 (a)	141316
Dense Drill - 4.0 (a)	141215

AMP | AMP-S

Application: Surgical
AMP-S and AMP
Ø3.75 or 4.0 (a)

Content:	
Tray	Code
Compact	131138

Implant Drivers	
Implant Insertion Driver - Medium	131139
Implant Insertion Driver - Medium-Manual	131146

Drills	
First Drill + 2.2	141138
Conical Drill - 2.8	141146
Conical Drill - 3.4 (a)	141148
Conical Drill - 3.8 (a)	141314
Dense Drill - 3.75 (a)	141316
Dense Drill - 4.0 (a)	141215

MT


Application: Surgical
MT Ø3.75

Content:	
Tray	Code
Compact	131138

Implant Drivers	
Implant Insertion Driver - Medium	131139
Implant Insertion Driver - Medium-Manual	131146

Drills	
First Drill + 2.2	141138
Conical Drill - 2.8	141146
Conical Drill - 3.4	141148
Dense Drill - 3.75	141316

Ø3.75 | Ø4.0





PROSTHETIC KIT

Prosthetic Kit


Application: All Interfaces and Components

Content:

Tray	Code
Compact	131138

Instruments

Torque Wrench - Surgical and Prosthetic	131128
Driver Adapter - Prosthetic Squared	131129
Hexagonal Driver - Short - 1.2	131010
Hexagonal Driver - Medium - 1.2	131011
Hexagonal Driver - Long - 1.2	131012
Conical Abutment Driver - Short	131016
Conical Abutment Driver - Medium	131017
Retriever - MT	131131



Includes Torque Wrench






Developed By P-I Brånemark

INSTRUMENT GUIDANCE-002-ENG-2019-03-13

www.pibranemark.com

 SIC invent Sweden AB
Stora Åvägen 21, 436 34
Askim - Sweden

 A COMPANY OF
SIC invent Group